

Class:

Roll No:


Experiment No:2  
To Determine the Wavelength of He-Ne LASER  
using a Diffraction Grating.

Date Of Preparation:

Date of Submission:

Signature of Teacher:

**Diagram:**



SCHEMATIC DIG.FOR He-Ne LASER DIFFRACTION  
THROUGH GRATING

## EXPERIMENT NO. 2

**AIM:** To determine the wavelength of the He-Ne LASER using a diffraction grating.

**APPARATUS:** LASER light source (He-Ne laser), diffraction grating and a meter scale

**THEORY:** LASER stands for Light Amplification by Stimulated Emission of Radiation. In He-Ne LASER the active medium is a mixture of He and Ne, and it is a four-level laser. A grating is an arrangement of a number of identical parallel slits. When a beam of LASER light is incident normally on the surface of a grating it gets diffracted to form a diffraction pattern made of zero order, first order, second order and higher orders depending on the grating element. The angle of diffraction ( $\Theta$ ) or each order (m) with respect to the direction of non-diffracted beam, the grating element (d) and the wavelength of incident light ( $\lambda$ ) are related by the equation known as GRATING EQUATION, viz.  $d \sin\Theta = m\lambda$

**OBSERVATION TABLE:** Grating Element (d) = 2.54/15000 cm

Sr. No.	Distance of grating from the screen D(cm)	Order of diffraction (m)	Distance between two maxima of same order (2X)	X (cm)	$\tan\Theta = X/D$
1					
2					
3					

**CALCULATIONS:**

Sr. No.	Angle of Diffraction ( $\Theta$ )	Sin ( $\Theta$ )	Wavelength $\lambda$ ( $\text{\AA}$ )	Mean $\lambda(\text{\AA})$
1				
2				
3				

**RESULTS:**

Calculated wavelength of He-Ne LASER = \_\_\_\_\_ ( $\text{\AA}$ )

Standard wavelength of He- Ne LASER = 6328 ( $\text{\AA}$ )

Percentage Error = \_\_\_\_\_