

Class:

Roll No:

Experiment No:1


To Determine the Wavelength of Sodium Light  
using a Diffraction Grating and Spectrometer.

Date Of Preparation:

Date of Submission:

Signature of Teacher:

Diagram



Schematic diag. for diffraction of Na light by  
Grating and Spectrometer

## EXPERIMENT NO. 1

**AIM:** To determine the wavelength of sodium light using a diffraction grating and spectrometer

**APPARATUS :** Diffraction grating, spectrometer, Sodium lamp, reading lens

**THEORY:** A plane diffraction grating is an arrangement of a large number of identical equidistant parallel slits. When a beam of monochromatic light is incident normally on the grating surface, it gets diffracted through various angles to give a diffraction pattern made of zero order, first order, second order and higher orders depending on the grating element  $d$ , (distance between two consecutive slits) and the wavelength of light being used. The grating equation relates the grating element ( $d$ ), the wavelength of incident light ( $\lambda$ ) and the angle of diffraction ( $\Theta$ ) for a given order ( $m$ ) and is given by

$$d \sin \Theta = m \lambda$$

**OBSERVATION TABLE:** Grating Element ( $d$ ) = (2.54/15000) cm

Order of maxima ( $m$ )	Spectrometer reading (L.H.S.) a (degrees) (MSR + VSR)	Spectrometer reading (R.H.S.) b (degrees) (MSR + VSR)	(a~b)= $2\Theta$	Angle of diffraction ( $\Theta$ )	Calculated wavelength ( $\lambda$ ) (Å)	Average wavelength ( $\lambda$ ) (Å)
0						
1						
2						

**CALCULATIONS:**             $d \sin \Theta = m \lambda$     (m-1)

$$d \sin \Theta = m \lambda \quad (\text{m-2})$$

**RESULTS :** The Average calculated value of wavelength = \_\_\_\_\_ (Å)

The standard value of wavelength = \_\_\_\_\_ (Å)

Percentage Error = \_\_\_\_\_ %